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4. CLINICAL PARTICULARS

4.1 Therapeutic indications:

For implantation. As a replacement for lost and insufficient tissue in the area

- Oral and maxillofacial surgery: Augmentation of the alveolar ridge in the maxilla, sinus lift, interpositional defects in Le Fort I osteotomy
- Reconstruction of the anterior cruciate ligament: Filling of tibia and femur tunnels
- Fractures of humerus/tibia: Proximal humerus fracture, depressed tibial plateau fractures
- Spinal surgery: Anterior cervical discectomy and fusion
- Hip arthroplasty: Septic/Revision Hip Arthroplasty

4.2 Posology and method of administration:

For grafting in accordance with the clinical indication. The type of application and localisation depends on the given surgical situation and is determined by the attending physician. The contents and outer surface of the inner bag are sterile. Prior to being used, human cortical substance, freeze-dried, DIZG, is to be rehydrated at room temperature in an appropriate sterile, physiological medium for at least 30 minutes (e.g. in isotonic infusion solution).

4.3 Contraindications:

Grafting in necrotic operation sites is contra-indicated. Due to an expected drop in the settling rate, accurate diagnosis of grafting in undersupplied or infected operation sites is imperative.

4.4 Special warnings and precautions for use:

For single use only.

4.5 Interactions with other medicinal products and other forms of interaction:

None known.

4.6 Pregnancy and lactation:

There are no sufficient data on the use of human corticalis, freeze-dried, DIZG during pregnancy or lactation.

4.7 Effects on ability to drive and operate machines:

Not applicable.

4.8 Undesirable effects:

None known.

The following incidence rates are used to assess side effects: very common (> 10 %); common (> 1 % - < 10 %); uncommon (> 0.1 % - < 1 %); rare (> 0.01 % - < 0.1 %); very rare (< 0.01 % or unknown).

Going on case history and physical examination results, each donor is carefully checked for any reasons for exclusion. Extensive serological tests for infection and tests for any -viruses are also performed using nucleic acid amplification techniques (NAT). The results of the donor's laboratory tests required for the procurement of the tissues (HIV1/2-Ab, HBC-Ab (IgG+IgM) (if positive, then anti-HBs > 10 I.E./l), HBsAg, HCV-Ab, TPHA, HTLV I/II-Ab and HBV-NAT, HCV-NAT, HIV-NAT (PCR) are negative. These tests exceed the requirements of EU Directive 2006/17. Despite these extensive tests and the use of a validated chemical cold sterilisation process, the transmission of infectious diseases through previously unknown or unidentifiable pathogens cannot be ruled out with complete certainty.

4.9 Overdose:

Not applicable.

5. PHARMACOLOGICAL PROPERTIES

Pharmacotheapeutic group:

Graft of human origin.

5.1 Pharmacodynamic and 5.2 Pharmacokinetic properties:

Owing to the primarily physical mechanism of action thereof, studies investigating the pharmacodynamics and pharmacokinetics of this therapeutic agent in trial subjects or patients are not available. Positive results in respect of the effectiveness and tolerability of the preparation in 288 patients have been substantiated in the „Study report on the clinical effectiveness and tolerability of allogenic, avital tissue implants sterilised in peracetic acid/ethanol“. This report is based on the statistically analysed and evaluated experience gained over a period of four years in the main hospitals using the preparation. It details excellent to good results in 98 % of all cases in respect of the use of human cortical substance as a spacer and restorative agent. Operating surgeons rated the settling as very good to good in 96 % of all cases. Positive results (4 patients) were achieved in the bridging-over of -defects in perioperative fractures among patients with status post implantation of titanium hip endoprotheses.

5.3 Preclinical safety data:

Given the decades of use in humans, experimental pre-clinical animal studies have not been conducted.

6. PHARMACEUTICAL PARTICULARS

6.1 List of excipients:

Not applicable.

6.2 Incompatibilities:

None known.

6.3 Shelf life:

5 years in the unopened container.

6.4 Special precautions for storage:

Store from + 4 °C to + 30 °C.

6.5 Nature and contents of container:

The contents and the external surface of the internal packaging are sterile.

Femur diaphysen: 1 piece

complete, 50 mm, 70 mm, 100 mm, 150 mm, 250 mm longitudinally halved

Femur diaphysis ring: 1 piece

10 mm, 12 mm, 14 mm, 16 mm

Femur chip: 1 piece

10 x 20 x 1 mm, 15 x 20 x 1 mm, 15 x 20 x 2 mm, 4 x 4 x 50 mm,

10 – 15 x 50 mm, 10 – 15 x 100 mm, 10 – 15 x 200 mm, 3-5 x 15 mm, 7-9 x 21 mm

Femur wedge: 1 piece

20 x 30 mm, 20 x 100 mm, 10 – 15 x 30,10 – 15 x 100 mm

Tibia diaphysen: 1 piece

complete, 50 mm, 70 mm, 100 mm, 150 mm, 250 mm longitudinally halved

Tibia chip: 1 piece

10 – 15 x 50 mm, 10 – 15 x 100 mm, 10 – 15 x 200 mm, 3-5 x 15 mm, 7-9 x 21 mm

Tibia wedge: 1 piece

10 – 15 x 30 mm, 10 – 15 x 100 mm

Fibula diaphysen: 1 piece

complete, 10 mm, 50 mm, 70 mm, 100 mm, 150 mm

Humerus diaphysis: 1 piece complete

Ulna diaphysis: 1 piece complete

Radius diaphysis: 1 piece complete

Costa: 1 piece complete, corpus

Granules:

Particle size 0.5 – 2.0 mm: 3,2 cm³, 6,4 cm³, 9,6 cm³

Particle size 0.5 – 1.25 mm: 3,2 cm³, 6,4 cm³, 9,6 cm³

Particle size 0.5 – 1.0 mm: 0,5 cm³, 1 cm³, 3 cm³, 5 cm³

Particle size 0,212 – 0,85 mm:

0,5 cm³, 1 cm³, 1,5 cm³, 3 cm³, 3,2 cm³, 5 cm³, 6,4 cm³, 9,6 cm³, 10 cm³, 15 cm³, 30 cm³

Particle size 1 – 2 mm: 0,5 cm³, 1 cm³, 1,5 cm³, 3 cm³, 5 cm³

Not all pack sizes may be marketed.

6.6 Special precautions for disposal and other handling:

The content of one pack is intended for use in one patient only. Any remaining tissue is to be discarded (disposal in accordance with AVV disposal key AS 180102). Upon being opened, the graft must be used within 2 hours.

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10. DATE OF REVISION OF THE TEXT

28.11.2025

11. GENERAL CLASSIFICATION FOR SUPPLY

Medicinal product subject to medical prescription.

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