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Please read carefully

1. NAME OF THE MEDICINAL PRODUCT:
Human cartilage tissue, cryopreserved, DIZG

2. COMPOSITION:
Proportions of human cartilage tissue.

3. PHARMACEUTICAL FORMS:
Cryopreserved graft pieces of various dimensions.

4. CLINICAL DATA

4.1 Therapeutic indications:
For implantation. For filling and mechanical stabilisation of tissue defects and as a placeholder for lost or insufficient tissue or for plastic surgery in the following specialist disciplines: General and plastic surgery, Thoracic, trauma and paediatric surgery, ENT medicine, orthopaedics, oral and maxillofacial surgery.

4.2 Dosage, method and duration of administration:
For transplantation according to the clinical indication. The graft must be thawed at room temperature for approximately 30 minutes in a suitable sterile physiological medium (e.g. an isotonic infusion solution) before use.

4.3 Contraindications:
Use in necrotic host tissue is contraindicated. For host sites that are poorly perfused or infected, the indication must be carefully restricted, as healing is less successful.

4.4 Special warnings and precautions for use:
For single use only.

4.5 Interactions with other medicinal products and other forms of interaction:
None known.

4.6 Pregnancy and breastfeeding:
There is insufficient information available on use during pregnancy or breastfeeding.

4.7 Influence on ability to drive and operate machinery:
Not applicable.

4.8 Undesirable effects:
None known.

The following frequencies are used to assess adverse reactions: Very common (>10%); Common (>1% - <10%); Uncommon (>0.1% - <1%); Rare (>0.01% - <0.1%); Very rare (<0.01% or unknown).

Each donor is examined for exclusion criteria on the basis of medical history and physical examination. This is followed by serological screening for infections, as well as viral testing using nucleic acid amplification techniques (NAT). The results of the laboratory tests required for the procurement of the donor's tissues for: HIV1/2 antibodies, HBc antibodies (IgG+IgM) (if positive, then Anti-HBs >10 IU/l), HBsAg, HCV antibodies, TPHA, HTLV I/II antibodies, and HBV-NAT, HCV-NAT, HIV-NAT (PCR) are negative. Despite these comprehensive tests, which go beyond the requirements of EU Directive 2006/17, and the use of a validated chemical cold sterilisation method, complete certainty that infectious diseases cannot be transmitted via previously unknown or undetectable pathogens cannot be guaranteed.

Reporting of suspected adverse reactions
Reporting suspected adverse effects after the authorisation of the medicinal product is important. It ensures the continuous monitoring of the medicinal product's benefit-risk profile. Healthcare professionals are asked to report any suspected adverse reactions to the Federal Institute for Vaccines and Biomedicines, Paul-Ehrlich-Institut, Paul-Ehrlich-Straße 51-59, 63225 Langen, Germany, phone: +49 6 10 37 70, fax: +49 61 03 77 12 34, website: www.pei.de. Patients should be advised to inform their doctor or healthcare provider if they notice any adverse reactions. This also applies to adverse reactions that are not listed in this Package Leaflet and the Summary of Product Characteristics. Patients may also report adverse reactions directly to the Paul Ehrlich Institute. By reporting adverse reactions, patients help increase knowledge about the safety of this medicinal product.

4.9 Overdose:
Not applicable.

5. PHARMACOLOGICAL PROPERTIES
Pharmacotherapeutic group:
Graft of human origin.

5.1 Pharmacodynamic and

5.2 Pharmacokinetic properties:
Due to the primarily physical mechanisms of action, there are no studies on pharmacodynamics or pharmacokinetics in patients or trial subjects.

5.3 Preclinical safety data:
As the material has been used in humans for decades, no preclinical animal studies have been carried out.

6. PHARMACEUTICAL PARTICULARS

6.1 List of excipients:
Not applicable.

6.2 Incompatibilities:
None known.

6.3 Shelf life:
2 years in unopened container.

6.4 Special precautions for storage:
Do not store above -40°C.

6.5 Type and contents of the container:
The contents and the outer surface of the inner packaging are sterile.

1 piece each
Rib cartilage 50 mm
Rib cartilage 100 mm
Articular cartilage

Not all pack sizes may be marketed.

6.6 Special precautions for disposal and other handling:
The contents of one pack are intended for use in a single patient only. Any remaining tissue residues must be discarded (disposal according to AVV waste code AS 180102). Once the container has been opened, the graft must be used within 2 hours.

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8. MARKETING AUTHORISATION NUMBER:
PEI.H.03360.01.1

9. DATE OF FIRST AUTHORISATION / RENEWAL OF THE AUTHORISATION:
23 December 2005 / 20 July 2010

10. DATE OF REVISION OF THE TEXT:
28 November 2025

11. GENERAL CLASSIFICATION FOR SUPPLY:
Prescription only

Note:
Please use the form on the DIZG homepage to report serious incidents or serious adverse reactions in accordance with Sections 8 and 9 TPG-GewV (Transplant Act-Tissue Ordinance).

Documentation:
Patient file labels are enclosed with each graft for documentation purposes under the TPG (Transplant Act).

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