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Package leaflet and Summary of Product Characteristics

Please read carefully

1. NAME OF THE MEDICINAL PRODUCT:

a) Fiberfill, lyophilised, DIZG

b) **Substance group:** Human tissue preparations, allogeneic.

2. THERAPEUTIC INDICATIONS

Congenital (e.g., cysts), traumatic or iatrogenic bone defects and fracture healing disorders (e.g., pseudarthrosis).

3. INFORMATION ON USE

a) Contraindications

Unstable spinal fractures; unstable endogenous bone structure; dental application; severe vascular disorders. Use in necrotic host sites is contraindicated. This medicinal product should be used only when strictly indicated in poorly perfused or infected host sites due to the poorer healing rate.

b) Precautions for use

For single use only.

Excessive rehydration of a graft can lead to loss of practicability or usability. Loss of practicability or usability may be indicative of excessive rehydration.

If too much rehydration solution is filled into the container, the excess liquid can be drained immediately after fixing the graft in place.

In the event of excessive rehydration: Allow the liquid to set, then drain the excess and manually homogenise. There is a risk of a fat embolism and/or embolism formation if the product enters a blood vessel.

Do not administer the product at high pressure and/or into a blood vessel!

c) Interaction with other medicinal products or other drugs that may potentially impact the effect of the medicinal product

None known.

d) Use in special patient populations: Pregnancy and lactation

There are no adequate data from use during pregnancy and breastfeeding.

e) Effects on ability to drive and use machines

Not applicable.

f) Special warnings and precautions for use

For single use only in one patient! Grafts must not be resterilised! Upon opening the container, the graft must be used within 2 hours.

4. INSTRUCTIONS FOR CORRECT USE

a) Posology in adults

For grafting in accordance with the clinical indication. The amount of Fiberfill to be used is dependent on the given surgical situation and is determined by the treating physician.

Prior to use, the graft is to be rehydrated with an appropriate physiological liquid, for example, physiological NaCl solution, which may, where indicated (e.g., septic bone defects), also contain an antibiotic. Glucose must not be used for reconstitution. For reconstitution, the required amount is poured into the recess of the graft after opening the screw-on container.

Pack size	Liquid to be added
5 cc	4.5 mL
10 cc	9 mL

The volumes indicated above correspond to when the recess is filled to capacity. After adding the liquid, the graft should be mixed in the vessel and brought to homogeneity. It is advised to wait until the transplant has absorbed the liquid and to apply it without excessive pressure.

The graft can be used 1 hour after adding the liquid.

b) Method of administration in adults

For implantation. Fiberfill is designed to fill defects when replacing lost or insufficient bone tissue. The method of administration and the site where Fiberfill is to be used is dependent on the given surgical situation and is determined by the treating physician.

c) Frequency of administration

The contents of one pack are intended exclusively for use in one patient. The amount of Fiberfill to be used and the frequency of administration are determined by the treating physician.

d) Duration of treatment

For single, permanent use.

e) Overdose: Symptoms, emergency procedures, antidotes

Not applicable.

5. UNDESIRABLE EFFECTS/ADVERSE REACTIONS

The following frequencies are used for the evaluation of undesirable effects: Very common (>10 %); Common (>1 % to <10 %); Uncommon (>0.1 % to <1 %); Rare (>0.01 % to <0.1 %); Very rare (<0.01 % or not known). Each donor will be checked for the presence of exclusion criteria on the basis of their medical history and physical examination. This will be followed by serological testing for infections, as well as testing for viral infections using the nucleic acid amplification technique (NAT). The results of the laboratory tests that are necessary to test the tissue harvested from the donor for: HIV-1/2-Ab, HAV-IgM-Ab or HAV-NAT, HBe-Ab (IgG+IgM) (if positive, then anti-HBs > 10 IU/L), HBsAg, HCV-Ab, TPHA, HTLV I/II-Ab and HBV-NAT, HCV-NAT, HIV-NAT (PCR) are negative. Despite these extensive tests, which exceed the requirements of EU Directive 2006/17, and the use of a validated chemical cold sterilisation process, the transmission of infectious diseases due to previously unknown or undetectable pathogens cannot be ruled out with complete certainty.

Pseudarthrosis, refracture, recurrence of the cyst, wound infection, haematoma, osteomyelitis, wound healing disorders, fat embolism / embolism, dysaesthesia, pain. Further surgical procedures may be required.

Reporting of suspected adverse reactions

Reporting suspected adverse reactions/side effects after authorisation of the medicinal product is important. It allows continued monitoring of the benefit/risk balance of the medicinal product. Healthcare professionals are asked to report any suspected undesirable effects/adverse reactions via the Federal Institute for Vaccines and Biomedicines, Paul-Ehrlich-Institut, Paul-Ehrlich-Straße 51 – 59, 63225 Langen, phone: +49 6 10 37 70, fax: +49 61 03 77 12 34, website: www.pei.de.

Patients must be instructed to tell their physician or the nurse if they notice any undesirable effects/adverse reactions. This also includes any undesirable effects/adverse reactions not listed in this package leaflet and Summary of Product Characteristics. Patients can also report undesirable effects/adverse reactions directly via the Paul-Ehrlich-Institut. By reporting undesirable effects/adverse reactions, patients contribute to providing more information on the safety of this medicinal product.

6. PHARMACOLOGICAL PROPERTIES

a) Pharmacodynamic,

b) Pharmacokinetic properties and

c) Toxicology: Based on the primarily physical mechanisms of action, no pharmacodynamic or pharmacokinetic studies in patients or trial subjects and no toxicology studies have been performed.

7. PHARMACEUTICAL PARTICULARS AND OTHER INFORMATION

a) Precautions for storage and shelf life

At +4 °C to + 30 °C the shelf life is 3 years in the unopened container.

b) Composition of the finished medicinal product (by type and quantity)

Fiberfill, lyophilised, DIZG, comprises a mixture of 25 % cancellous granules and 75 % demineralised, human cortical fibres.

Fiberfill contains 55.2 mg or 110.4 mg phosphate buffer in a mixture of sodium chloride, disodium hydrogen phosphate, potassium dihydrogen phosphate.

c) Date of dispensing: see consignment note.

d) Pharmaceutical form: Lyophilised mixture of cancellous granules and demineralised, cortical fibres. Graft of human origin

Nature and contents of container by weight, volume or number of units:

The content, the closed container, the inner aluminium bag as well as the inside of the outer bag are sterile.

Per unit: 5 cm³, 10 cm³.

Not all pack sizes may be marketed.

e) Details of the marketing authorisation holder

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[German Institute for Cell and Tissue Replacement]

Köpenicker Straße 325, 12555 Berlin;

Phone: +49(0)30-6576-3050; fax: +49(0)30-6576-3055;

Email: dizg@dizg.de

f) Details of the manufacturer who authorised marketing of the finished medicinal product

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[German Institute for Cell and Tissue Replacement]

Köpenicker Straße 325, 12555 Berlin;

Phone: +49(0)30-6576-3050; fax: +49(0)30-6576-3055;

Email: dizg@dizg.de

g) Approval number: PEI.G.12020.01.1

h) Date of first authorisation: 08/02/2021

i) Medicinal product status: Prescription-only medicine.

8. OTHER INFORMATION

Precautions for handling and disposal:

Any remaining tissue must be discarded.

Disposal in accordance with German waste disposal laws [AVV Abfallschlüssel AS 180102)].

9. DATE OF REVISION OF THE TEXT

18/10/2022

Transplantation consignment note:

Every tissue graft is supplied with a transplantation consignment note. Any adverse drugs reactions (particularly infections) that occur in the immediate postoperative period, as well as their circumstances, are to be recorded on the enclosed report form. The transplantation consignment note serves to record and evaluate the clinical results as well as any.

DIZG

German Institute for Cell and Tissue Replacement

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