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1. NAME OF THE MEDICINAL PRODUCT
Human cancellous bone, freeze-dried, DZIG

2. COMPOSITION
Human cancellous bone tissue.

3. PHARMACEUTICAL FORMS
Cubes, chips, granules, cylinders.

4. CLINICAL PARTICULARS
4.1 Therapeutic Indications:
For implantation as tissue replacement for lost or insufficient tissue and for filling and stabilizing bony defects. It is used in the field of general, neuro and accident surgery as well as in orthopaedics. Further applications in other operative disciplines are also conceivable.

4.2 Posology and method of administration:
For grafting in accordance with the clinical indication. The contents and outer surface of the inner bag are sterile. Prior to being used, human cancellous bone, freeze-dried, DZIG, is to be rehydrated at room temperature in an appropriate aseptically, physiological medium for at least 30 minutes (e.g. in isotonic infusion solution).

4.3 Contraindications:
Grafting in necrotic operation sites is contra-indicated. Due to an expected drop in the settling rate, accurate diagnosis of grafting in undersupplied or infected operation sites is imperative.

4.4 Special warnings and precautions for use:
For single use only.

4.5 Interactions with other medicinal products and other forms of interaction:
None known.

4.6 Pregnancy and lactation:
There are no sufficient data on the use of human cancellous, freeze-dried, DZIG during pregnancy or lactation.

4.7 Effects on ability to drive and operate machines:
Not applicable.

4.8 Undesirable effects:
None known.

The following incidence rates are used to assess side effects: very common (>10%); common (>1% – <10%); uncommon (>0.1% – <1%); rare (>0.01% – <0.1%); very rare (<0.01% or unknown).

The following incidence rates are used to assess side effects: very common (>10%); common (>1% – <10%); uncommon (>0.1% – <1%); rare (>0.01% – <0.1%); very rare (<0.01% or unknown).

The following incidence rates are used to assess side effects: very common (>10%); common (>1% – <10%); uncommon (>0.1% – <1%); rare (>0.01% – <0.1%); very rare (<0.01% or unknown).

The following incidence rates are used to assess side effects: very common (>10%); common (>1% – <10%); uncommon (>0.1% – <1%); rare (>0.01% – <0.1%); very rare (<0.01% or unknown).

4.9 Overdose:
Not applicable.

5. PHARMACOLOGICAL PROPERTIES
Pharmacotherapeutic group:
Graft of human origin.

5.1 Pharmacodynamic and 5.2 Pharmacokinetic properties:
Owing to the primarily physical mechanism of action thereof, studies investigating the pharmacodynamics and pharmacokinetics of this therapeutic agent in trial subjects or patients are not available. Positive results in respect of the effectiveness and tolerability of the preparation in 1205 patients have been substantiated in the „Study report on the clinical effectiveness and tolerability of allogenic, avial tissue implants sterilised in peracetic acid/ethanol“. This report is based on the statistically analysed and evaluated experience gained over a period of four years in the main hospitals using the preparation. It details excellent to good results in 90% of all cases in respect of the use of human cancellous as a spacer and restorative agent. Operating surgeons rated the settling as very good to good in 82% of all cases. Positive results (4 patients) were achieved in the bridging-over of defects in perioperative fractures among patients with at least post implantation of titanium hip endoprosthesis.

5.3 Preclinical safety data:
Given the decades of use in humans, experimental pre-clinical animal studies have not been conducted.

6. PHARMACEUTICAL PARTICULARS

6.1 List of excipients:
Not applicable.

6.2 Incompatibilities:
None known.

6.3 Shelf life:
5 years in the unopened container.

6.4 Special precautions for storage:
Store from +4 °C to +30 °C.

6.5 Nature and contents of container:
The contents and the external surface of the internal packaging are sterile.

cancellous cubes/cylinders:

| | |
|--|--|
| 1 piece | 1 cm ³ , 3 cm ³ , 9 cm ³ , 12 cm ³ , 15 cm ³ |
| 3 pieces | 1 cm ³ , 3 cm ³ |
| 3 cm ³ cancellous chips 1 – 2 mm | 3 cm ³ , 5 cm ³ , 15 cm ³ , 30 cm ³ |
| cancellous chips 2 – 4 mm | 3 cm ³ , 5 cm ³ , 15 cm ³ , 20 cm ³ , 30 cm ³ |
| Corticocancellous chips <1 cm: 30 x 20 mm | 1 cm ³ , 3 cm ³ , 5 cm ³ , 10 cm ³ , 15 cm ³ , 30 cm ³ , 45 cm ³ , 60 cm ³ |
| Corticocancellous chips 1, shape 15 x 20 mm, 10 mm thick | 1 piece |
| Corticocancellous wedges: 10 x 10 – 15 mm, 4 – 9 mm thick | length 25 – 30 mm, height and width 15 mm |
| Granules: | |
| Particle size 0.212 – 0.85 mm: | 0.5 cm ³ , 1 cm ³ , 3 cm ³ , 5 cm ³ |
| 0.5 mm, 1.0 mm, 1.5 mm: | 0.5 cm ³ , 1.0 cm ³ , 1.5 cm ³ |
| Particle size 1 – 2 mm: | 1 piece |
| cancellous cylinder: | 15 x 15 mm (height), 12 x 30 mm (height), 15 x 15 mm (height), 20 x 15 mm (height) |
| 10 x 30 mm (height), 12 x 30 mm (height), 15 x 15 mm (height), 20 x 15 mm (height) | 1 piece and 2 pieces |
| Caput femoris: | 1 piece and 2 pieces |
| endoprosthesis, halved, without cartilage: | 1 piece and 2 pieces |
| endoprosthesis, quartered, without cartilage: | 1 piece |
| halved, without cartilage: | 1 piece |
| quartered, without cartilage: | 1 piece |
| Wedge: | 1 piece, 22.5°, 45° |
| Os ilium, bilateral: | 1 piece |
| Os ilium, unilateral: | 1 piece |
| 10 – 19 x 40 mm, 10 – 19 x 60 mm, 20 x 20 mm, 20 x 30 mm, 20 x 40 mm, 20 x 60 mm, 20 x 80 mm, 30 x 20 mm, 40 x 20 mm, 60 x 20 mm | Not all pack sizes may be marketed. |
| 10 – 19 x 40 mm, 10 – 19 x 60 mm, 20 x 20 mm, 20 x 30 mm, 20 x 40 mm, 20 x 60 mm, 30 x 20 mm, 40 x 20 mm, 60 x 20 mm | |
| Not all pack sizes may be marketed. | |

6.6 Special precautions for disposal and other handling:
The content of one pack is intended for use in one patient only. Any remaining tissue is to be discarded (disposal in accordance with AVV disposal key AS 180120). Upon being opened, the graft must be used within 2 hours.

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10. DATE OF REVISION OF THE TEXT
18.04.2019

11. GENERAL CLASSIFICATION FOR SUPPLY
Medicinal product subject to medical prescription.

Adverse outcomes attributable to the tissue must be promptly reported to your local representative.

Documentation: patient file labels are enclosed with each graft for documentation purposes.

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