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**3. PHARMACEUTICAL FORM**

Diaphyses, chips, wedges, granules.

**4. CLINICAL PARTICULARS**
**4.1 Therapeutic indications:**

For implantation to bridge over bone defects. Further applications in other operative disciplines are also conceivable.

**4.2 Posology and method of administration:**

For grafting in accordance with the clinical indication. The contents and outer surface of the inner bag are sterile. Prior to being used, human cortical substance, freeze-dried, DIZG, is to be rehydrated at room temperature in an appropriate sterile, physiological medium for at least 30 minutes (e.g. in isotonic infusion solution).

**4.3 Contraindications:**

Grafting in necrotic operation sites is contra-indicated. Due to an expected drop in the settling rate, accurate diagnosis of grafting in undersupplied or infected operation sites is imperative.

**4.4 Special warnings and precautions for use:**

For single use only.

**4.5 Interactions with other medicinal products and other forms of interaction:**

None known.

**4.6 Pregnancy and lactation:**

There are no sufficient data on the use of human corticalis, freeze-dried, DIZG during pregnancy or lactation.

**4.7 Effects on ability to drive and operate machines:**

Not applicable.

**4.8 Undesirable effects:**

None known.

The following incidence rates are used to assess side effects: very common (>10%); common (>1% – <10%); uncommon (>0.1% – <1%); rare (>0.01% – <0.1%); very rare (<0.01% or unknown).

Going on case history and physical examination results, each donor is carefully checked for any reasons for exclusion. Extensive serological tests for infection and tests for any viruses are also performed using nucleic acid amplification techniques (NAT). The results of the donor's laboratory tests required for the procurement of the tissues (HIV1/2-Ab, HAV-IgM-Ab or HAV-NAT, HBc-Ab (IgG+IgM) (if positive, then anti-HBs >10 I.E./l), HBsAg, HCV-Ab, TPHA, HTLV I/II-Ab and HBV-NAT, HCV-NAT, HIV-NAT (PCR)) are negative. These tests exceed the requirements of EU Directive 2006/17. Despite these extensive tests and the use of a validated chemical cold sterilisation process, the transmission of infectious diseases through previously unknown or unidentifiable pathogens cannot be ruled out with complete certainty.

**4.9 Overdose:**

Not applicable.

**5. PHARMACOLOGICAL PROPERTIES**
**Pharmacotherapeutic group:**

Graft of human origin.

**5.1 Pharmacodynamic and 5.2 Pharmacokinetic properties:**

Owing to the primarily physical mechanism of action thereof, studies investigating the pharmacodynamics and pharmacokinetics of this therapeutic agent in trial subjects or patients are not available. Positive results in respect of the effectiveness and tolerability of the preparation in 288 patients have been substantiated in the "Study report on the clinical effectiveness and tolerability of allogenic, avital tissue implants sterilised in peracetic acid/ethanol". This report is based on the statistically analysed and evaluated experience gained over a period of four years in the main hospitals using the preparation. It details excellent to good results in 98% of all cases in respect of the use of human cortical substance as a spacer and restorative agent. Operating surgeons rated the settling as very good to good in 96% of all cases. Positive results (4 patients) were achieved in the bridging-over of defects in peri-operative fractures among patients with status post implantation of titanium hip endoprotheses.

**5.3 Preclinical safety data:**

Given the decades of use in humans, experimental pre-clinical animal studies have not been conducted.

**6. PHARMACEUTICAL PARTICULARS**
**6.1 List of excipients:**

Not applicable.

**6.2 Incompatibilities:**

None known.

**6.3 Shelf life:**

5 years in the unopened container.

**6.4 Special precautions for storage:**

Store from +4 °C to +30 °C.

**6.5 Nature and contents of container:**

The contents and the external surface of the internal packaging are sterile.

Femur diaphyses: 1 piece

complete, 50 mm, 70 mm, 100 mm, 150 mm,

250 mm longitudinally halved

Femur diaphysis ring: 1 piece

10 mm, 12 mm, 14 mm, 16 mm

Femur chip: 1 piece

10 x 20 x 1 mm, 15 x 20 x 1 mm, 15 x 20 x 2 mm, 4 x 4 x 50 mm,

10 – 15 x 50 mm, 10 – 15 x 100 mm, 10 – 15 x 200 mm, 3 – 5 x 15 mm, 7 – 9 x 21 mm

Femur wedge: 1 piece

20 x 30 mm, 20 x 100 mm, 10 – 15 x 30,10 – 15 x 100 mm

Tibia diaphyses: 1 piece

complete, 50 mm, 70 mm, 100 mm, 150 mm

250 mm longitudinally halved

Tibia chip: 1 piece

10 – 15 x 50 mm, 10 – 15 x 100 mm, 10 – 15 x 200 mm, 3 – 5 x 15 mm, 7 – 9 x 21 mm

Tibia wedge: 1 piece

10 – 15 x 30 mm, 10 – 15 x 100 mm

Fibula diaphyses: 1 piece

complete, 10 mm, 50 mm, 70 mm, 100 mm, 150 mm

Humerus diaphysis: 1 piece complete

Ulna diaphysis: 1 piece complete

Radius diaphysis: 1 piece complete

Costa: 1 piece complete, corpus

Granules:

Particle size 0.5 – 2.0 mm: 3,2 cm<sup>3</sup>, 6,4 cm<sup>3</sup>, 9,6 cm<sup>3</sup>Particle size 0.5 – 1.25 mm: 3,2 cm<sup>3</sup>, 6,4 cm<sup>3</sup>, 9,6 cm<sup>3</sup>Particle size 0.5 – 1.0 mm: 0,5 cm<sup>3</sup>, 1 cm<sup>3</sup>, 3 cm<sup>3</sup>, 5 cm<sup>3</sup>Particle size 0.212 – 0.85 mm: 0,5 cm<sup>3</sup>, 1 cm<sup>3</sup>, 1,5 cm<sup>3</sup>, 3 cm<sup>3</sup>, 3,2 cm<sup>3</sup>, 5 cm<sup>3</sup>, 6,4 cm<sup>3</sup>, 9,6 cm<sup>3</sup>, 10 cm<sup>3</sup>, 15 cm<sup>3</sup>, 30 cm<sup>3</sup>Particle size 1 – 2 mm: 0,5 cm<sup>3</sup>, 1 cm<sup>3</sup>, 1,5 cm<sup>3</sup>, 3 cm<sup>3</sup>, 5 cm<sup>3</sup>

Not all pack sizes may be marketed.

**6.6 Special precautions for disposal and other handling:**

The content of one pack is intended for use in one patient only. Any remaining tissue is to be discarded (disposal in accordance with AVV disposal key AS 180102). Upon being opened, the graft must be used within 2 hours.

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**10. DATE OF REVISION OF THE TEXT**

24.01.2024

**11. GENERAL CLASSIFICATION FOR SUPPLY**

Medicinal product subject to medical prescription.

Adverse outcomes attributable to the tissue must be promptly reported to your local representative.

**Documentation:** patient file labels are enclosed with each graft for documentation purposes.

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Customer:

Leaflet:

Item number / Print number (customer material number):

Format open:

Format closed:

Code:

Font size / line spacing:

Colour:

Prepress editor and signature:

## Digital proofs are not colour consistent!

We will not be held responsible for any errors which are not marked on this proof. Author's corrections (subsequent changes to text, colours, format etc.) are invoiced separately and lead to a change of the expected delivery date.

We need your approval to start the actual print production.

**Please note that we need your approval 15 working days before the delivery date.**

Please sign, date and return

**Approved**

**Resubmit after correction**

Date, signature

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Edelmann Burscheid GmbH