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Human fascia connective tissue.

## 3. PHARMACEUTICAL FORMS

Freeze-dried grafts of various sizes.

## 4. CLINICAL PARTICULARS

### 4.1 Therapeutic indications:

For implantation. Positive results have been achieved in medial collateral ligament stabilisation during knee operations. Further applications in other operative disciplines are also conceivable, e.g. as tissue replacement or spacer for lost or insufficient tissue or for plastic surgery.

### 4.2 Posology and method of administration:

For grafting in accordance with the clinical indication. Prior to being used, human fascia, freeze-dried, DIZG, is to be rehydrated at room temperature in an appropriate sterile, physiological medium for at least 30 minutes (e.g. in isotonic infusion solution).

### 4.3 Contraindications:

Grafting in necrotic operation sites is contra-indicated. Due to an expected drop in the settling rate, accurate diagnosis of grafting in undersupplied or infected operation sites is imperative.

### 4.4 Special warnings and precautions for use:

For single use only.

### 4.5 Interactions with other medicinal products and other forms of interaction:

None known.

### 4.6 Pregnancy and lactation:

There are no sufficient data on the use of human fascia, freeze-dried, DIZG during pregnancy or lactation.

### 4.7 Effects on ability to drive and operate machines:

Not applicable.

### 4.8 Undesirable effects:

None known.

The following incidence rates are used to assess side effects: very common (>10%); common (>1%–<10%); uncommon (>0.1%–<1%); rare (>0.01%–<0.1%); very rare (<0.01% or unknown).

Going on case history and physical examination results, each donor is carefully checked for any reasons for exclusion. Extensive serological tests for infection and tests for any viruses are also performed using nucleic acid amplification techniques (NAT). The results of the donor's laboratory tests required for the procurement of the tissues (HIV1/2-Ab, HBC-Ab (IgG+IgM) (if positive, then anti-HBs > 10 I.E./l), HBsAg, HCV-Ab, TPHA, HTLV I/II-Ab and HBV-NAT, HCV-NAT, HIV-NAT (PCR)) are negative. These tests exceed the requirements of EU Directive 2006/17. Despite these extensive tests and the use of a validated chemical cold sterilisation process, the transmission of infectious diseases through previously unknown or unidentifiable pathogens cannot be ruled out with complete certainty.

### 4.9 Overdose:

Not applicable.

## 5. PHARMACOLOGICAL PROPERTIES

### Pharmacotherapeutic group:

Graft of human origin.

### 5.1 Pharmacodynamic and 5.2 Pharmacokinetic properties:

Owing to the primarily physical mechanism of action thereof, studies investigating the pharmacodynamics and pharmacokinetics of this therapeutic agent in trial subjects or patients are not available. Positive results in

respect of the effectiveness and tolerability of the preparation in 50 patients have been substantiated in the „Study report on the clinical effectiveness and tolerability of allogenic, avital tissue implants sterilised in peracetic acid/ethanol“. This report is based on the statistically analysed and evaluated experience gained over a period of four years in the main hospitals using the preparation. It details positive results in 100% of all cases in respect of the use of human fascia lata for collateral ligament/capsular instability of the medial knee joint.

### 5.3 Preclinical safety data:

Given the decades of use in humans, experimental pre-clinical animal studies have not been conducted.

## 6. PHARMACEUTICAL PARTICULARS

### 6.1 List of excipients:

Not applicable.

### 6.2 Incompatibilities:

None known.

### 6.3 Shelf life:

3 years in the unopened container.

### 6.4 Special precautions for storage:

Store from +4 °C to +30 °C.

### 6.5 Nature and contents of container:

The contents and the external surface of the internal packaging are sterile.

*Fascia lata:* 1 piece each

20 x 70 mm, 20 x 100 mm, 20 x 200 mm, 20 x >250 mm,

40 x 50 mm, 40 x 80 mm, 40 x 100 mm, 40 x 200 mm,

15 x 20 mm, 20 x 30 mm, 30 x 40 mm.

Not all pack sizes may be marketed.

### 6.6 Special precautions for disposal and other handling:

The content of one pack is intended for use in one patient only. Any remaining tissue is to be discarded (disposal in accordance with AVV disposal key AS 180102). Upon being opened, the graft must be used within 2 hours.

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## 10. DATE OF REVISION OF THE TEXT

28.11.2025

## 11. GENERAL CLASSIFICATION FOR SUPPLY

Medicinal product subject to medical prescription.

Adverse outcomes attributable to the tissue must be promptly reported to your local representative.

**Documentation:** patient file labels are enclosed with each graft for documentation purposes.

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